



What you should know about

Juzo Compression Garments



Freedom in Motion

For nearly 100 years, our family philosophy has been that medical garments should not inhibit a person's lifestyle.

Frequently Asked Questions:

When should I wear a garment?

First thing in the morning and remove at the end of the day.

When should I purchase a new compression garment?

Every 6 months or as therapy progresses and new symptoms develop.

Can I use lotions & ointments? What about Latex allergies?

Juzo® compression garments are **LATEX FREE** using high quality Lycra® covered with textile fibers that offer wearing comfort, air permeability, and are largely resistant to sun, ointments, lotions and perspiration.

Where can I find material content, size and style?

A tag sewn inside provides details of the garment.

What's inside a Juzo compression package?

Care & application brochure highlights use and related information. Open toe stockings include a donning device (open toe slippie) and instructions.

What does graduated compression mean?

Compression begins at the ankle and decreases up the leg. For arm sleeves, the compression begins at the wrist and decreases.

Care and Maintenance:

How to wash your garment?

Washing helps restore the garment back to it's original form.

Hand or machine wash warm & dry daily.

Use the lowest permanent press setting.

Garment bag is recommended.

Separate dark colors from other clothing.

Use Juzo® washing detergent or another mild soap.

DO NOT USE: bleach, chlorine, fabric softeners, or other laundry additives.

*Always use Juzo latex free donning gloves when applying compression garments.
Rough fingernails, jagged rings etc. can damage fabric.*

For more information, go to:

www.juzousa.com

Patient Therapy Help Line: 888-255-1300 • support@juzousa.com