

SOLARIS
FREEDOM TO LIVE

*New Glove
& Gauntlet
Styles*

EXOTM
BY SOLARIS

*Off-the-Shelf Gradient
Compression Garments*

Toomac  **Solutions**
DIVISION OF TOOMAC HOLDINGS LTD

32C Poland Rd, Glenfield, Auckland 0627
www.toomac.co.nz

Ph: 09 443 5347
Fax: 09 443 5345



Exo garments provide a full range of ready-to-wear solutions for your patients' unique needs. From ExoStrong ready-to-wear flat-knit garments for maximum containment to ExoFusion for focused compression over the foot, there's an Exo for everyone.



Why Choose Exo?



Exo garments are easy to care for - machine wash and dry with a normal load of laundry.



ExoStrong is the complete line of ready-to-wear flat-knit garments and only costs a fraction of custom daywear.



Choose the style and structure that best accommodates your patients' unique needs.



Based on feedback from our customers, we've updated the target compression on ExoStrong garments. It still provides your customers with 20 - 30 mmHg, but now it will be easier for them to don.



ExoSoft arm sleeves were the highest-rated circular-knit products by therapists and DMEs in single blind trials*.

*Visit solarismed.com/exo_comparison to see the full results.

ExoStrong™ Upper Extremity Garments



ExoStrong Arm Sleeve 54-221

A ready-to-wear arm sleeve designed for maximum support, comfort, and durability from the wrist to the axilla.



Designed with a functional bend at the elbow for improved comfort and range of motion.



Flat-knit structure provides maximum containment while remaining flexible.



Durable material helps resist snags and runs.



Made with fibers that combine a soft, tactile feel with moisture wicking and temperature control.

Options		
Colors		
Beige	Black	
Proximal Options		
	Knit Top	Silicone Top
Compression		
20 - 30 mmHg		

ExoStrong Glove / Gauntlet 54-203 / 54-211

A comfortable, off-the-shelf glove / gauntlet that provides maximum containment.



Flat, reinforced edges help reduce irritation and rolling.



Tapered wrist design works in tandem with UE-CG garments to provide full arm coverage.



Gusseted web spaces provide flexibility while maintaining containment.



Made with fibers that combine a soft, tactile feel with moisture wicking and temperature control.

Options	
Colors	
Beige	Black
Compression	
20 - 30 mmHg	
Available Styles	
Full Finger Glove DIPs to Wrist	
Half Finger Glove PIPs to Wrist	
Quarter Finger Glove MCPs to Wrist	
Gauntlet MCPs to Wrist	


Pictured:
Quarter Finger
MCPs to Wrist

For more information, visit solarismed.com

Care Instructions

Washing

- Machine wash with mild detergent in warm water.
- DO NOT use bleach, fabric softeners, dryer sheets, or other additives.



Drying

- Machine dry on permanent press.
- DO NOT iron.



Indications

- Edema

Contraindications

- Arterial insufficiency or degeneration
- Untreated congestive heart failure
- Systemic infection
- Impaired sensation
- Cognitive impairment
- Allergies or sensitivities to materials

For more information, visit solarismed.com