

Arion donning and doffing aids

Product benefits

Easy to Use



- Minimize friction between compression garment and skin
- Reduce physical strain for both compression wearers and caregivers

Safe to Apply



- Reduce risk of skin damage ; can be especially beneficial in case of high skin sensitivity or leg ulcers
- May contribute to a longer garment lifetime by preventing damage from fingers and nails
- Help maintaining the original compression longer and avoiding garments from overstretching and distortion

Simple Solutions



- For many types of compression garments: stockings, arm sleeves, and gloves
- Designed for different shapes and sizes of limbs
- Support the compliance to compression therapy

Care instructions



Product overview

Download our complete Arion application aids catalogue:



Compression
therapy made
easier



Arion donning and doffing aids
for medical compression garments

for Kids

Toomac Solutions
DIVISION OF TOOMAC HOLDINGS LTD

Distributed by:
Toomac Solutions
32D Poland Rd, Glenfield
Auckland 0627
Freephone 0508 443 5347
Email: enquiries@toomac.co.nz
Wed: www.toomacsolutions.co.nz

ARION
SLIDE SOLUTIONS

Arion Slide Solutions B.V.
Rutherford 11 . 6422 RE Heerlen . NL
Tel.: +31 45 209 4014
E-Mail: slideorders@arion-group.com
www.slide-solutions.com



Follow us on: [@arionslidesolutions_en](https://www.instagram.com/arionslidesolutions_en)

T552z-011-EN (From 01/REV 09-2023)

ARION
SLIDE SOLUTIONS

Arion donning aids


for children



Putting on medical compression garments for children is a challenge for both the children and their parents or carers, and it can be very painful for the little patients. However, in order to achieve a successful therapy, this daily recurring action is absolutely necessary. The Arion donning aids for children have therefore been specially developed to simplify the donning of medical compression garments such as gloves, arm sleeves, stockings or tights in the case of burn injuries, scalds or scar healing.

Wearing compression garments inhibits the development of scars and helps prevent further scarring. In addition, the compression garment helps the scars to heal faster and better.

Our application aids are made of a particularly smooth material, which means that there is little friction on the child's skin when putting on the garment. Performing this daily task therefore takes less effort for the caregiver and is less stressful for the little patients than without using a donning aid.




ARION Easy-Slide Kids



On

Donning aid for medical compression arm sleeves and stockings with an open toe



Art. Nr.	Size	Length
E01700	Extra Small	29,5 cm / 11,5"
E01701	Small	39,5 cm / 15,5"




Video instructions

ARION Magnide Kids



On

Donning aid for medical compression stockings with a closed toe



Art. Nr.	Size	Shoe Size
E01712	Small	≤ 36 / ≤ 3,5



Video instructions




ARION Dactyna Kids



On

Donning aid for medical compression gloves with open fingertips



Video instructions

Art. Nr.	Size	Finger circumference	Finger lengths
E01770	Extra Small	≤ 5,5 cm / 2,2"	6 cm / 2,4"



Easy-Grip Gloves



On



Off



On



On

Anti-slip gloves as an additional donning and doffing aid in combination with the other Arion devices

These gloves provide the user an additional grip and help put on compression stockings, arm sleeves and compression gloves more gently. They prevent from pulling on the compression garment and protect from damages caused by jewelry or fingernails. This contributes to preserve the material and to extend the lifespan of the compression garment.

The Arion anti-slip gloves contain latex.



Art. Nr.	Size	Hand size EU	Hand size UK
E08030	Small	< 7,5 cm	< 3"
E08010	Medium	7,5 - 8,5 cm	3 - 3,5"
E08020	Large	8,5 - 9,5 cm	3,5 - 4"
E08040	Extra large	> 9,5 cm	> 4"

