

Arion Application Gloves



Necessary accessory for donning and doffing medical compression garments in combination with Arion Application Aids.

Specialty designed latex-free gloves provide an optimal grip on the fabric and help medical stockings to slide over the foot and up the leg.

Now you can avoid pulling hard at the top of your stockings, thus preserving the original medical compression and ensuring a longer life of the garment. Furthermore, Arion Application Gloves protect the fabric of stockings from being damaged by nails or jewelry.



Size	Hand size EU (cm)	Hand size UK (inch)
Small	< 7,5 cm	< 3 inch
Medium	7,5 - 8,5 cm	3 - 3,5 inch
Large	8,5 - 9,5 cm	3,5 - 4 inch
Extra large	> 9,5 cm	> 4 inch



**We make your
life easier!**

Toomac Solutions
DIVISION OF TOOMAC HOLDINGS LTD

Toomac Solutions


32D Poland Rd
Wairau Valley
Auckland
Phone 0800 443 5347
Email: sales@toomac.co.nz

TSS19-022-2

Arion Export B.V.
Rutherford 11
6422 RE Heerlen, The Netherlands
Tel.: +31 (0)45 - 209 40 14
E-mail: exp@arion-group.com
www.arion-group.com

**Easy on -
Easy off**

ARION
simply smile



Arion Application Aids
for your medical
compression garments

Easy to Use

- Minimizes friction between compression garment and skin thanks to the product's unique sliding properties
- Reduces daily stress associated with putting on and taking off medical compression
- Allows the users of medical compression to apply their garments by themselves


sim-slide®

easy, safe and simple

The optimal multi-functional application aid for donning and doffing open toe medical compression stockings and tights in one product.

Sim-Slide® makes an application of medical stockings easier, thus reducing daily stress associated with donning and doffing compression garments.

A Study shows that, when using Arion application aids, medical compression wearers achieve greater satisfaction* with successfully putting on and taking off their compression garments. This contributes to their positive user experience, thus improving therapy compliance and optimizing therapy results.



Size	Shoe size
Small	< 34 - 36
Medium	37 - 38,5
Large	39 - 41
Extra large	≥ 42



Safe to Apply

- Lowers the risk of skin damage; this can be especially beneficial for care recipients with high skin sensitivity or leg ulcers
- Helps to maintain clinical benefits of original medical compression by avoiding garment overstretching/distortion
- Ensures longer life of compression garments without damage from fingers or nails

magnide® 2in1

easy, safe and simple

The new and improved application aid for putting on and taking off closed toe medical compression stockings and tights in one product.

Thanks to the product's unique sliding properties Magnide® 2in1 offers a smooth transition of the garment onto the patient's skin. It makes the application of compression garment safe and easy, minimizing its wear and tear.

Magnide® 2in1 has flexible sizing to fit different sizes of limbs. This is particularly important for patients diagnosed with lower limb oedema (swelling).

Magnide® 2in1 is a simple solution for healthcare professionals to recommend and for care recipients to use in combination with medical compression therapy.



Size	Shoe size
Medium	≤ 38
Large	≥ 39 f/m 43
Extra large	≥ 44
XL also suitable for extra wide feet with CA upto 28 cm	



Simple Solution

- Offers application aids for all types of medical compression: stockings, sleeves, gloves
- Designed for different shapes and sizes of limbs
- Empowers concordance of medical compression wearers to apply their garments regularly

easy-slide arm

easy, safe and simple

The user friendly application aid for putting on medical compression arm sleeves

The Arion Application Aid Easy-Slide Arm makes putting on compression arm sleeves quick and easy.

A straightforward donning method in combination with the application aid requires very little physical effort on the part of the compression wearer or his caregiver.

Simply put the aid over your arm, and slide your medical compression sleeve up the limb over the application aid.

Easy-Slide Arm helps compression wearers enjoy the convenience of putting on compression arm sleeves independently and regularly.



Size	Circumference of widest part of the arm
Medium	37 - 38,5
Large	39 - 41

