

Easy-Grip Gloves




Anti-slip gloves to put on and take off compression garments gently and with a good grip, as a complement to the Arion application aids

These gloves have been specially developed to enable you to put on compression stockings over the foot and the leg by providing an optimal grip. With the Arion anti-slip gloves, you avoid pulling hard at the border of the stockings, preserving their material and contributing to a longer lifespan for the garment. In addition, the gloves protect the fabric of your stockings from being damaged by jewelry or fingernails.

The Arion anti-slip gloves contain latex.

Care instructions:



Art. Nr.	Size	Hand size EU	Hand size US
Eo8030	Small	< 7,5 cm	< 3"
Eo8010	Medium	7,5 - 8,5 cm	3 - 3,5"
Eo8020	Large	8,5 - 9,5 cm	3,5 - 4"
Eo8040	Extra large	> 9,5 cm	> 4"

Anti-Slip Mat



Anti-slip mat to facilitate putting on and taking off compression stockings, as a complement to the Arion application aids

- Creates a soft anti-slip surface
- Makes it easier to position the stocking correctly on the heel
- Lightweight and easy to handle

Art. Nr: Eo8003
Size: 37,5 x 19,5 x 0,5 cm /
14.8" x 7.7" x 0,2"

Care instructions:

Clean with gentle soap and water,
air dry



Arion donning and doffing aids

PRODUCT BENEFITS

Easy to Use



- Minimize friction between compression garment and skin
- Reduce physical strain for both compression wearers and caregivers
- Promote independence of compression wearers

Safe to Apply



- Reduce risk of skin damage ; can be especially beneficial in case of high skin sensitivity or leg ulcers
- May contribute to a longer garment lifetime by preventing damage from fingers and nails
- Help maintaining the original compression longer and avoiding garments from overstretching and distortion

Simple Solutions



- For many types of compression garments: stockings, tights, arm sleeves, and gloves
- Designed for different shapes and sizes of limbs*
- Support the compliance to compression therapy

*Made-to-measure product possible upon request: ask your reseller or Arion!

Toomac Solutions
DIVISION OF TOOMAC HOLDINGS LTD

Distributed by:
Toomac Solutions
32D Poland Rd, Glenfield
Auckland 0627
Freephone 0508 443 5347
Email: enquiries@toomac.co.nz
Web: www.toomacsolutions.co.nz

ARION
SLIDE SOLUTIONS

Arion Slide Solutions B.V.
Rutherford 11 . 6422 RE Heerlen . NL
Tel.: +31 45 209 4014
E-Mail: slideorders@arion-group.com
www.slide-solutions.com



Follow us on: @arionslidesolutions_en

TSS23-04 (Toomac) (REV 05-2023)

Compression therapy

made easier



Arion donning and doffing aids
for medical compression garments

ARION
SLIDE SOLUTIONS



Two solutions in one product: donning and doffing aids for compression stockings and tights

Thanks to their unique sliding properties, Magnide® 2in1 and Sim-Slide® reduce friction between the leg and the medical compression stocking. As a result, the stocking can be easily put on without pulling.


Our application aids are designed for different shapes and sizes of limbs. This is especially important for compression wearers with oedema in the lower legs.

Magnide® 2in1 and Sim-Slide® are simple solutions for specialists and compression wearers, making medical compression therapy more convenient.

For stockings with a closed toe



On + Off



Art. Nr.	Size	Shoe Size (UK)	
		Woman	Man
E01540	Medium	≤ 38	≤ 6
E01541	Large	39 - 43	6,5 - 9
E01542	X-Large	≥ 44	≥ 9,5

XL also suitable for extra wide feet (circumference up to 28 cm / 11")



Video instructions:






On

Extension aid to put on medical compression stockings and tights in combination with Magnide® (2in1), Easy-Slide and Sim-Slide®

Arion SlideX covers the leg and reduces friction and pain from above the ankle to the knee and even higher.*

Arion SlideX offers an easy solution for compression wearers and caregivers for putting on compression stockings and tights in case of:

- open wounds with a plaster or bandage**
- leg ulcers
- a larger than average lower leg (e.g. lipoedema/ lymphoedema)
- high stockings or pantyhose*



Art. Nr.	Size	Circumference of the leg***
E01520	Medium	40 cm / 15.7"
E01521	Large	60 cm / 23.6"
E01522	Extra large	100 cm / 39.4"

*Sizes L and XL are longer and cover the leg above the knee: these are the right sizes to help put on thigh-high stockings or tights.

**Open wounds should be covered with an appropriate dressing before using the product.

***Of the widest part of the leg



Video instructions:




On

The user-friendly donning aid for medical compression arm sleeves

The donning aid Arion Easy-Slide Arm makes donning compression arm sleeves as easy as child's play.


Only little effort is required on the part of the user or his helper. Simply pull the application aid over the arm and let the medical compression arm sleeve slide over it. Easy-Slide Arm supports the independence of the compression wearer, as the compression garment can be put on without assistance, which contributes to therapy compliance.



Art. Nr.	Size	Arm length*	Max. product circumference**
E01505	Medium	59 - 67 cm / 23" - 26.5"	45.6 cm / 18"
E01506	Large	> 67 cm / > 26.5"	57.6 cm / 22.5"

* Arm length from the tip of the middle finger to the armpit.

** For a full-length distribution of the Easy-Slide Arm over your arm, we recommend that the circumference of the widest part of your arm is min. 1 cm/0.4" smaller than the max. product circumference. If this is not possible with the current sizes available, please ask your reseller or Arion for a made-to-measure product.



Video instructions: