Accessories



Juzo ACS Light Hip Band

Art. J711


The Juzo ACS Light Hip Band can be used for fixing the thigh wrap, if required.



Juzo ACS Light Midline Band

Art. J712


The Juzo ACS Light Midline Band has been specially developed for the compression of the torso (e.g. thorax, pelvis, abdomen and/or genital area). It can be cut to size and is easy to put on, thanks to the general-purpose velcro fasteners.



Juzo ACS Light Velcro Fastener Sheet

Art. J713

With the Juzo ACS Light Velcro Fastener Sheet, additional velcro fasteners can be individually cut to size



Juzo ACS Light Pads

Art. J714

Juzo ACS Light Pads can be used to provide extra wearing comfort with additional cushioning and/or to provide increased pressure in specific areas.



Juzo ACS Light Anti-slip Band Art. J715

The Juzo ACS Light Anti-slip Band offers an additional fastening option and can be shortened as required. It can be worn under the midline band and hip band as well as under the thigh wrap to provide secure fastening.




Juzo ACS Light Extender Art. J716

The Juzo Light Extender provides increased compression for the extremities (e.g. the ankle, back of the hand or instep).









Juzo ACS Light

Adjustable compression system

Juzo ACS Light short-stretch compression segments are used for additional treatment during decongestive therapy. The products are easy to put on, thanks to the incorporated lining (fixing aid). They can be adapted to the individual patient using adjustable fasteners – even if circumferences change during the therapy. This uncomplicated application increases patient's compliance and supports self-treatment. With Juzo ACS Light, therapeutic compression up to 40 mmHg can be achieved.

Features

- + high working pressure and low resting pressure
- + compression levels up to 40 mmHg
- + particularly wear-resistant material
- + can be worn day and night
- + adaptable whilst wearing for ideal compression
- + easy to apply yourself even with limited mobility or strength
- + easy to put on by the patient, thanks to the sewn-in fixing aid
- + various optional accessories
- + Juzo ACS light segments can also be combined with Juzo SoftCompress and Juzo Classic Seamless products

Indications

- · lymphoedema
- · lipoedema
- · phleboedema
- post-traumatic oedema
- · postoperative oedema

Colours





BI

Upper extremity



Lower extremity



Juzo ACS Light Thigh Wrap Art. J704LI

- + with sewn-in silicone border
- + with incorporated lining (fixing aid)
- + loose extender strap
- + available with extender in left and right versions



- + with sewn-in silicone border
- with incorporated lining (fixing aid)
- + with knee comfort hole
- + can be worn both left and right



Juzo ACS Light Calf Wrap Art. J702LI

- + with sewn-in silicone border
- + with incorporated lining (fixing aid)
- + available in left and right versions



Juzo ACS Light Foot Wrap Art. J701LI

- + with incorporated malleo pads
- + loose extender strap
- + can be worn both left and right

Highlights



Velcro fasteners

Velcro fasteners make it easy to apply the product and to adjust it in case of fluctuations in body circumference without taking the product down.



Secure hold

Thanks to the incorporated lining (fixing aid), the segments can be put on in next to no time and provide a perfect fit. The sewn-in silicone border prevents slipping.



Wearing comfort

hing A separate pad included with the be arm and hand wrap offers additioovi- nal comfort.



Anatomical fit

The compression wraps are perfectly adapted to individual body contours. The wraps for the lower leg, thigh and arm are available for both the left and right side of the body.