



<section-header><section-header>

Freedom in Motion

Quality standards Verified quality



GÜTEZEICHEN

MEDIZINISCHE KOMPRESSIONS

STRÜMPFE

RAL

Oeko-Tex® Standard 100

"Oeko-Tex Standard 100" is a globally-standardised and independent testing and certification system that applies to all stages of processing, from raw materials through to the intermediate and final product. The more intensive the skin contact, the stricter the requirements that must be met.

All Juzo compression fabrics, both the circular knitted and flat knitted varieties, have been awarded the "Oeko-Tex Standard 100" label for textiles. This enables the company Juzo to guarantee that all components that go into a product are completely harmless to humans and the environment. The overall requirements of this standard far exceed those of existing national legislation.

Quality seal RAL


All medical compression stockings from Juzo comply with the relevant RAL standard. The test specifications contain a comprehensive set of binding guidelines on the manufacture of medical compression stockings and sleeves, for example:

- material characteristics (harmless to humans)
- periodic quality inspections, for example by the Hohenstein Research Institute
- rubber-elastic properties in lengthwise and crosswise directions
- compression classes 1 4, for sleeves compression classes 1 3
- manufacture of custom-made and standard sized stockings



ISO certification

Juzo implemented its general quality management system back in 1998, and was one of the first companies in this industry to be certified to EN ISO 9001 and EN 46001. Juzo is currently certified to EN ISO 9001 and EN ISO 13485, which superseded EN 46001.



Medical compression garments The 4 pressure levels

Medical compression stockings from Juzo guarantee correct pressure values and a well-defined pressure gradient. The pressure decreases continuously from the ankle section upwards towards the heart, and also towards the toes. We can provide a wide variety of garments, from below-knee stockings to pantyhose, from gloves to face masks, depending on the therapeutic application.

Compression fabrics are categorised into 4 compression classes (CCL) based on the level of compression in the ankle or wrist area according to the quality seal RAL. Possible indications can be, for example:

18 - 21 mmHg (CCL 1) - mild compression

- $\cdot\,\,$ minor obstructive disorders accompanied by feelings of tiredness or heaviness in the legs
- · isolated varicose vein formation
- · preventative measure during pregnancy
- 23 32 mmHg (CCL 2) moderate compression
- varicose veins
- · leg swelling during pregnancy
- after varicose vein sclerosis or operations
- · for the long-term treatment of healed leg ulcers
- · propensity for oedema formation exists

34 – 46 mmHg (CCL 3) – strong compression

- · varicose veins accompanied by damage to the deep leg veins
- · after healing of venous crural ulcers
- severe swelling tendency of the legs
- chronic venous insufficiency
- · severe propensity for oedema formation
- > 49 mmHg (KKL 4) sehr starke Kompression
- severe lymphatic obstructions
- · lymphoedema with fibrosis
- · lymphoedema stage 3

Pregnancy The body in an exeptional state

During pregnancy, a woman's body has great demands placed on it. Up until the birth of the child, the blood volume in the mother increases by 30 - 40 %. This is very taxing on the veins. Furthermore, the hormones released in the body soften the connective tissue to prepare the body for the changes to come. This also has the effect of making the vein walls more elastic and they can lose their tension.

The growing uterus also presses on the inferior vena cava. This is the largest vein in the body, which carries blood from the lower limbs, parts of the pelvis, and from the stomach organs. Pressure on this vein impedes the return flow of blood, which increases the risk of spider veins or varicose veins. Pregnant women are therefore advised to counteract this risk through movement and regular exercise. Many pregnant women complain of "swollen legs", especially in summer. Don't be afraid to take a break now and then and elevate your legs. Your body will thank you for it. Compression stockings also provide relaxation and relief for heavy and swollen legs. They can be worn for the entire duration of the pregnancy, and you can start wearing them from early pregnancy onwards. This will help you enjoy your pregnancy to the fullest!

· Party and the state of the state of the





Juzo Inspiration Elegance and uniqueness without limits

Colours			Compression classes	1 2 3 4
Sugar	Ту	/pe	Model	Juzo Inspiration Art. 2701, 2702
Sesame	Le	engths in cm		
Almond	A	D	Below-knee stockings	
Cinnamon	SI	ormal 40 – 46 nort 34 – 39 ktra short 31 – 34	Top border Thin silicone border · 3.5 cm Thin balance silicone border · 3.5 cm	•
Cacao	A	F	Half-thigh stockings	
Poppy seed Blueberry	4		Silicone border · 5 cm Silicone border pattern · 5 cm Balance silicone border · 5 cm Balance silicone border pattern · 5 cm Lace silicone border · 5 cm	•
Black pepper	A	G	Thigh stockings	
Trend Colours ¹ For current Trend Colours see www.juzo.com/trend-colours Fashion Colours ¹	SI	ormal 72 – 83 nort 66 – 72 ktra short 58 – 66	Silicone border / wide · 5 cm Silicone border pattern / wide · 5 cm Balance silicone border / wide · 5 cm Balance silicone border pattern / wide · 5 cm Lace silicone border · 5 cm	•
For current Fashion Colours	A	г	Pantyhose	Hip size 80 – 120 cm
bip Dye Collection ¹	N Sł	side leg measurement ormal 72 – 83 nort 66 – 72 «tra short 58 – 66	Closed crotch	•
Available in colour Poppy seed and Blueberry as half-thigh	A	TU	Pregnancy pantyhose	Art. 2781, 2782
stockings (AF), thigh stockings (AG) and pantyhose (AT). For details see www.juzo.com/dip-dye Juzo Inspiration	ln N	regnancy body part side leg measurement ormal 72 – 83 nort 66 – 72	Closed crotch	•
Batik Collection ¹		ktra short 58 – 66		
All standard colours (except Sugar and Almond) and Trend Colours can be batiked in black		pen		•
or white www.juzo.com/batik-collection	CI	osed		•
Vision pattern		• Ava	ilable in standard sizes I — VI and custom-made • Onl	y available in custom-made

Vision pattern is available in standard sizes with closed toes

¹Production time 8 days, non-returnable



Page 6



Juzo Soft

Activity boosts your mood and your well-being

Colours		Compression classes	1 2 3 4
Sugar ¹	Type Lengths in cm	Model	Juzo Soft Art. 2001, 2002
Sesame	AD	Below-knee stockings	A.C. 2001, 2002
Almond	Normal 40 – 46 Short 34 – 39 Extra short 31 – 34	Top border Thin silicone border · 3.5 cm Thin balance silicone border · 3.5 cm	•
Cinnamon	AF	Half-thigh stockings	
Cacao		Top border Silicone border · 5 cm Silicone border pattern · 5 cm Balance silicone border · 5 cm Balance silicone border pattern · 5 cm Lace silicone border · 5 cm	•
Blueberry	AE with his attachment		
Black pepper	AF with hip attachment	Half-thigh stocking Stocking left, right or "to be worn as a pair"	•
Juzo Soft	AG	Thigh stockings	
Trend Colours ¹ For current Trend Colours see www.juzo.com/trend-colours Batik Collection ¹	Normal 72 – 83 Short 66 – 72 Extra short 58 – 66	Top border Silicone border / wide - 5 cm Silicone border pattern / wide - 5 cm Balance silicone border / wide - 5 cm Balance silicone border pattern / wide - 5 cm Lace silicone border - 5 cm	•
All standard colours (except	AG with hip attachment	Thigh stocking	
Sugar and Almond) and Trend Colours can be batiked in black or white		Stocking left, right or "to be worn as a pair"	•
www.juzo.com/batik-collection	AT	Pantyhose	
Production time 8 days, non-returnable	Hip size 80 – 140 cm Inside leg measurement Normal 72 – 83 Short 66 – 72 Extra short 58 – 66	Closed crotch Open crotch Fly	•
_	BT	Leggings	
X	Hip size 80 – 140 cm Inside leg measurement Normal 72 – 83 Short 66 – 72	Closed crotch	CCL 1, not in custom- made
	Toes		
	Open		•
L	Closed Short foot ³		•

• Available in standard sizes I – VI and custom-made • Only available in custom-made

Silicone borders			
Silicone border	Balance silicone border	Lace silicone border	
Silicone border pattern	Balance silicone border pattern		

Additional Products

Art.	Model	Art.	Model
9800	Juzo Adhesive Lotion	1701	Arion Easy-Slide Arm
i Dana	Special roll-on adhesive for comfortable wearing of compression stockings and compression sleeves Contents each: 60 ml		Donning aid for compression sleeves and thorax segments Smooth, two-layer material, very easy to use Sizes: M, L
9830	Juzo Vital Balsam 7	1712	Arion Easy-Slide Kids
C Juzo	Special skin care for compression therapy Cools, nourishes and revitalises irritated skin Contents each: 100 ml Contents each: 30 ml		Donning aid for compression sleeves, stockings and thorax segments Especially for children Smooth, two-layer material, very easy to use
9860	Juzo Lymph Lotion	4700	Sizes: XS, S
03.00	Skin care lotion with 4 % urea and 30 % lipids for daily use in the treatment of lymphoedema and lipoedema and also venous afflictions Contents each: 500 ml	1702	Arion Easy-Slide Leg Donning aid for compression stockings and pantyhose with open toes Smooth, two-layer material, very easy to use Sizes: S, M, L, XL
9000	Juzo Special Detergent	1705	Arion Sim-Slide®
osut 6	For medical compression garments, all types of elastic fabrics and textiles, and high quality fabrics made from wool and silk Highly concentrated and economical to use Contents each: 250 ml		Donning and doffing aid for compression stockings and pantyhose with open toes Smooth, two-layer material, simply fasten to the ankle using the elastic strap Sizes: S, M, L, XL
	(sufficient for 130 I wash water) Contents each: 30 mI	1703	Arion Magnide®
9301	(sufficient for 15 I wash water) Juzo Special gloves		Donning aid for compression stockings and pantyhose with closed toes
Ein	High-quality textile glove with rubberised palm for easy donning Very durable	1709	Smooth material, very easy to use due to innovative magnetic closure Sizes: M, L, XL Arion Magnide [®] 2in1
~~ ~	Sizes: S/5, M/7, L/9, XL/11	1705	Donning and doffing aid for compression
1890	Juzo Easy Fit Donning and doffing aid for compression stockings with open and closed toes Sizes: M, XL	1711	stockings and pantyhose with closed toes Smooth material, very easy to use due to innovative magnetic closure Sizes: M, L, XL Arion Easy-Off®
1900	Juzo Donning aid		Doffing aid for compression stockings and
	Donning aid for compression stockings and pantyhose with open toes	J	pantyhose with open and closed toes Sizes: S, M, L
	Pliable nylon fabric, smooth surface	1708	Arion Dactyna
1901	Juzo Slippie Leg Donning aid for compression stockings and pantyhose with open toes Made of strong Tyvek, excellent gliding effect		Donning aid for compression gloves with open fingers Smooth, two-layer material, very easy to use Sizes: XS, S, M, L
1902	Juzo Slippie Arm Donning aid for compression sleeves and thorax segments Tearproof material excellent gliding effect		

 $\textbf{Toomac Solutions} \cdot 32 \text{D} \text{ Poland } \text{RD} \cdot \text{Wairau Valley} \cdot \text{Auckland} \cdot \text{New Zealand} \cdot \text{juzo} @ \text{toomac.co.nz} \cdot \text{www.toomacsolutions.co.nz} \\ \textbf{Solutions} \cdot \textbf{Solution$

 $\textbf{Julius Zorn GmbH} \cdot Julius platz 1 \cdot 86551 Aichach \cdot Deutschland \cdot info@juzo.de \cdot www.juzo.de$