



VETERINARY KINESIOLOGY TAPE

PRODUCT LINE



WWW.VETKINTAPE.COM

MEDICAL TAPING



Medical Taping is recognised as a proven and internationally accepted method in physiotherapy clinics and hospitals all over the world. It provides positive, reliable results and is expected to remain one of the treatment methods used daily by many physiotherapists. THYSOL teamed up with veterinary professionals to develop a more robust tape, specifically for use on animals, giving the same results as the human variant, CureTape®. The development of VetkinTape® was possible thanks to our in-house production and R&D facilities.

VETKINTAPE®

VetkinTape® has been developed, in collaboration with veterinary professionals, exclusively for animal welfare, and is produced in our certified factory in South Korea. It is an elasticated cotton, latex-free tape giving good ventilation to the treatment site. VetkinTape® has a hypoallergenic acrylic adhesive layer, which offers 25% increased adhesiveness compared to existing tapes on the market. Unlike our human CureTape®, VetkinTape® is not suitable for human skin.



In practice we have seen that VetkinTape® can be used to assist in:

1. Injured muscles – the tape lifts the skin from the underlying tissues. This can reduce pain and, by promoting circulation, assists in the reduction of any swelling by enhancing blood flow to the injured muscle.
2. Biomechanical dysfunction – the tape provides support and stability without restricting the range of motion.
3. Postural dysfunction – the tape improves the proprioceptive feedback to the area thereby assisting in retraining the muscles to maintain correct posture.



Medical Taping can be used either as a sole or as a supplementary treatment in animal therapy throughout the rehabilitation process.

VETKINTAPE® IS USED FOR:

- Activate blood and lymph circulation.
- Treatment of biomechanical dysfunction.
- Pain redemption.
- Influencing fascia and muscles.
- Treatment of scars.
- Injury lesions and inflammations.
- (After) treatment of injuries to the musculoskeletal device.
- Improvement of muscle function.
- Treatment of postural dysfunction.

To ensure the optimal application of VetkinTape®, THYSOL offers canine and equine professional courses for veterinary professionals across the world. These courses have been developed in co-operation with equine and canine practitioners from a veterinary and physiotherapy background, and other professionals from the Animal Health industry. There is continuing close co-operation with several veterinary universities and education organizations.

VetkinTape® is available in 50 countries worldwide and is seen as the global leader in veterinary taping.



VETKINTAPE 3 CM

The 3 cm wide VetkinTape is widely used in treatment of smaller muscle groups, and is recommended for smaller animals like dogs.

Size 1 roll: 3 cm x 5 m
Packaging unit: 20 rolls per display box (10 boxes with 2 rolls)
Available colours: ● ● ● ●

VetkinTape 3 cm x 5 m Blue	Article no.: 2082
VetkinTape 3 cm x 5 m Red	Article no.: 2083
VetkinTape 3 cm x 5 m Black	Article no.: 2084
VetkinTape 3 cm x 5 m Orange	Article no.: 2085



VETKINTAPE 6 CM

The 6 cm wide VetkinTape is the most used size in clinics. Also commonly used for cutting lymph strips and other treatment specific tape shapes such as X and Y tape.

Size 1 roll: 6 cm x 5 m
Packaging unit: 10 rolls per display box
Available colours: ● ● ● ●

VetkinTape 6 cm x 5 m Blue	Article no.: 1082
VetkinTape 6 cm x 5 m Red	Article no.: 1083
VetkinTape 6 cm x 5 m Black	Article no.: 1084
VetkinTape 6 cm x 5 m Orange	Article no.: 1085




VETKINTAPE 10 CM

The 10 cm wide VetkinTape is suitable for broad muscle groups. Also commonly used for cutting lymph strips and other tape shapes such as X and Y tape.

Size 1 roll: 10 cm x 5 m
Packaging unit: 6 rollen per displaydoos
Available colours: ● ● ● ●

VetkinTape 10 cm x 5 m Blue	Article no.: 3082
VetkinTape 10 cm x 5 m Red	Article no.: 3083
VetkinTape 10 cm x 5 m Black	Article no.: 3084
VetkinTape 10 cm x 5 m Orange	Article no.: 3085





LITERATURE

This instruction book gives a powerful insight into the many and varied possibilities for using VetkinTape. It contains practical, clearly illustrated examples of more than 50 conditions which can be treated with a range of different taping techniques.

MTC book "Kinesiology Taping For Horses"
by Renate Ettl - ISBN / EAN: 978-90-825826-0-4.
Content: 207 pages with 276 full colour photos.

Article No. : 2081

CROSSLINQ®

CrossLinqs are synthetic, non-elastic, raster-shaped tapes with a latex-free adhesive layer. CrossLinqs are applicable on pain points, trigger points, scar tissue, acupuncture points and they affect the energy fields in the body (meridians).

CrossLinq Small: 20 sheets x 12 pieces (1.4 x 1.8 cm) Article no.: 1502
CrossLinq Medium: 20 sheets x 8 pieces (2.1 x 2.7 cm) Article no.: 1503
CrossLinq Large: 20 sheets x 4 pieces (4,4 x 5,2 cm) Article no.: 1504

SCISSORS

VetkinTape scissors; Soft Touching Article no.: 4020
VetkinTape scissors; High Quality Article no.: 4040

SPRAY

The Pre-Taping spray is used to clean and degrease the coat before applying the tape.

Clean Coat Pre-Taping spray, 200 ml. Article no.: 7010

LEATHER ACCESSORIES

To handle scissors and tape during the treatment, the following leather accessories have been developed. Multiple roller holders can be attached to the belt.

VetkinTape leather belt with leather roll and scissor holder Article no.: 4501
VetkinTape leather scissor holder Article no.: 4504
VetkinTape leather roll holder (for 1 roll; 6 cm or 2 x 3 cm) Article no.: 4505



THYSOL Group BV
Midzomerweg 8
7532 SW Enschede
The Netherlands
www.thysol.com



WWW.VETKINTAPE.COM