



THYSOL®



**Toomac**  **Solutions**  
DIVISION OF TOOMAC HOLDINGS LTD

PROFESSIONAL TAPE  
PRODUCT CATALOGUE

# ABOUT THYSOL® GROUP

THYSOL is the market leader throughout Europe, supplying its products to both the medical industry and patients alike. THYSOL GROUP is among the very few companies that has control over the production facilities. This ensures THYSOL has complete control during the manufacturing process guaranteeing high quality production. This also allows THYSOL to be extremely innovative in supplying a variety of new products.

These standards and the brands below are the foundation of the company;

**CureTape®**, a portfolio of medical level (human) kinesiology tapes.

**VetkinTape®**, a product line of veterinary kinesiology tape.

**FASCIQ®**, a range of fascia massage tools

Having conducted over 15,000+ courses worldwide to date, THYSOL continues to train medical professionals on how to use their products through the MTC® brand.

**MTC®** or Medical Training Concept is a series of courses that specialize in teaching health care professionals how to best use THYSOL's extensive range of health aids to the benefit of the patient via a series of detailed educational courses and workshops.

## Our Vision

At THYSOL, our aim is to enable both therapist and patient alike to gain easy access to affordable innovative treatment methods and high-quality products to prevent "immobility" and promote optimum health.

## Our Mission

As a leading tape company, we provide the facilities, products and trained personnel to deliver state-of-the-art medical devices and services that enable people and animals to live life to the fullest.

Our philosophy is that our products in combination with the right application should help, support and speed-up the recovery of the body.

# CONTENT

<b>CureTape®</b> .....	<b>4</b>
CureTape® Classic .....	6
CureTape® Big Boy Classic .....	6
CureTape® Giant Classic.....	7
CureTape® Punch .....	7
CureTape® Sports .....	8
CureTape® Giant Sports.....	8
CureTape® Big Boy Sports.....	9
CureTape® Art.....	9
CureTape® Beauty.....	10
CureTape® Gentle .....	10
CureTape® 3 Roll Cylinder .....	11
CureTape® Books.....	12
CureTape® Accessories .....	13
<b>CrossLinq®</b> .....	<b>14</b>
CrossLinq® Patches .....	14
CrossLinq® Book.....	14
<b>FASCIQ®</b> .....	<b>16</b>
FASCIQ® IASTM Tools.....	18
FASCIQ® Cupping Sets .....	20
FASCIQ® Tape.....	21
FASCIQ® Accessories .....	21

# ABOUT CURETAPE®

**CureTape®** elastic kinesiology tapes are skin-friendly and anti-allergic. They are latex free, certified with the TÜV quality mark and produced in the ISO-13485/9001 certified factory in South Korea.

The high quality elastic cotton is water resistant and is backed with a breathable, hypoallergenic acrylic adhesive layer. CureTape products are Medical Class 1 registered and fulfill all European (EUDAMED – MDR 2020), FDA and international medical quality requirements.

Due to the high level of quality and control, CureTape is safe to use and therefore the best choice for the medical professional to achieve the finest result possible within the treatment.

## Basic principles for using CureTape® For medical treatment

- It activates the body's natural recovery process.
- (Post)treatment of musculoskeletal injuries.
- Improvement of muscle function and skin structures.
- Full freedom of movement.
- It's a 24/7 treatment and allows tissue to recover optimally after trauma.

## Effects

The elasticity of CureTape in relation to the elasticity of the skin is used to create a lifting effect on the epidermis. The visual effect of lifting the skin with tape are called "convolutions".

CureTape stimulates the fascial structures through the skin and has an effects on all kind of systems in the body. It activates the body's own recovery mechanism, with optimal freedom of movement.

By using different taping techniques, a range of different effects can be achieved.

## Characteristics

- Water-permeable, while still being water resistant.
- Elasticity of the tape is similar to that of the human skin (130%-140%).
- Only stretchable in length.
- High-quality ventilating material.
- The tape can be worn for several days (3-5 days).
- Developed for professional medical use.
- Medical class 1 registered.
- Latex free with TÜV quality mark.
- Anti-allergic acrylic adhesive backing applied in wave-form.
- Skin friendly, no added substances.





# CureTape® Classic



**Materials:** Elastic cotton with a breathable, hypoallergenic acrylic adhesive.

**Dimensions:** Roll: 5 cm x 5 meters

**Displaybox:** 6 rolls

**Available colours:**

CureTape kinesiotape is the most used brand among professionals. The Classic variant is our bestseller: versatile, skin-friendly and meets the highest quality requirements. CureTape Classic is an elastic, breathable, cotton tape with a hypoallergenic acrylic adhesive layer. A bond time of 4 to 7 days can be expected for most tape applications.



**Classic Beige**  
8717624160103



**Classic Blue**  
8717624160134



**Classic Black**  
8717624160226



**Classic Pink**  
8717624160165



**Classic Green**  
8717624161216



**Classic Orange**  
8717624160196



**Classic Yellow**  
8717624161087



**Classic Red**  
8717624161254

### Additional Sizes



**7.5 cm x 5 m**  
**Classic Beige**  
8717624160257



**2.5 cm x 5 m**  
**Classic Beige**  
8717624160288



**1 cm x 5 m**  
**Classic Beige**  
8717624160929



**1 cm x 5 m**  
**Classic Blue**  
8717624160950



# CureTape® Big Boy Classic

ON REQUEST

**Materials:** Elastic cotton with a breathable, hypoallergenic acrylic adhesive.

**Dimensions:** Roll: 5 cm x 5 meters

**Displaybox:** 40 rolls (10 rolls per colour)

**Available colours:**

**Big Boy Classic**  
8717624163005



# CureTape® Giant Classic



**Materials:** Elastic cotton with a breathable, hypoallergenic acrylic adhesive.

**Dimensions:** Roll: 5 cm x 31 meters

**Displaybox:** 1 roll

**Available colours:**

The CureTape Giant is the largest roll in our range, measuring 5 cm x 31.5 m. CureTape Giant classic has the same composition as CureTape Classic with an elastic, breathable, cotton tape with a hypoallergenic acrylic adhesive layer. Most tape applications have a bond time of 4 to 7 days.



**Giant Classic Beige**  
8717624160318




**Giant Classic Blue**  
8717624160332



**Giant Classic Black**  
8717624160394



**Giant Classic Pink**  
8717624160356



**Giant Classic Green**  
8717624161308



**Giant Classic Orange**  
8717624160370



**Giant Classic Yellow**  
8717624161049



**Giant Classic Red**  
8717624161322

# CureTape® Punch



**Materials:** Elastic cotton with a breathable, hypoallergenic acrylic adhesive.

**Dimensions:** Roll: 5 cm x 5 meters

**Displaybox:** 6 rolls

**Available colours:**

Due to the unique hole pattern, this tape has a greater stretching capacity and the working direction is three-dimensional. The skin and fascia receive much more signals due to the pressure differences that arise around the hole pattern. The skin also has more opportunity to ventilate.



**Punch Beige**  
8717624160622



**Punch Blue**  
8717624160653



**Punch Black**  
8717624160745



**Punch Pink**  
8717624160684

# CureTape® Sports



**Materials:** Elastic viscose tape with a breathable, hypoallergenic acrylic layer with 25% more adhesive strength than the Classic range.

**Dimensions:** Roll: 5 cm x 5 meters

**Displaybox:** 6 rolls

**Available colours:**

CureTape Sports kinesio tape is an elastic tape made of viscose with 25% more adhesive strength. The perfect choice for sports-related complaints and injuries, but also for other applications where more adhesive strength is desired, such as tape applications on feet and hands. The sports variant is also often chosen for preventive use during sports. CureTape Sports, is made of viscose. Viscose is just as comfortable to wear as cotton, it has a smoother, more supple feel to it, it dries faster and is easy to apply.



**Sports Beige** 8717624164033  
**Sports Blue** 8717624162701  
**Sports Black** 8717624162671  
**Sports Pink** 8717624165320  
**Sports Lime** 8717624162640  
**Sports Orange** 8717624162169

# CureTape® Giant Sports



**Materials:** Elastic viscose tape with a breathable, hypoallergenic acrylic layer with 25% more adhesive strength than the Classic range.

**Dimensions:** Roll: 5 cm x 31 meters

**Displaybox:** 1 roll

**Available colours:** **ON REQUEST**

The CureTape Giant is the largest roll in our range, measuring 5 cm x 31.5 m. CureTape Giant sports has the same composition as CureTape sports with an elastic, breathable, viscose tape tape with a hypoallergenic acrylic adhesive layer with 25% more adhesive power. Most tape applications have a bond time of 4 to 7 days.



**Giant Sports Beige** 8717624165702  
**Giant Sports Blue** 8717624165719  
**Giant Sports Black** 8717624165726  
**Giant Sports Pink** 8717624165733  
**Giant Sports Lime** 8717624165740  
**Giant Sports Orange** 8717624165757

# CureTape® Big Boy Sports

**ON REQUEST**



**Materials:** Elastic viscose tape with a breathable, hypoallergenic acrylic layer with 25% more adhesive strength than the Classic range.

**Dimensions:** Roll: 5 cm x 5 meters

**Displaybox:** 40 rolls (10 rolls per colour)

**Available colours:**



**Big Boy Sports**  
8717624165788

# CureTape® Art




**Materials:** Elastic cotton with a breathable, hypoallergenic acrylic adhesive.

**Dimensions:** Roll: 5 cm x 5 meters

**Displaybox:** 6 rolls

**Available colours:**

The Art variant is a real eye-catcher with its modern, colourful prints; versatile, skin-friendly and meets the highest quality requirements. The various prints can be beautifully combined with the uni colors of CureTape Classic. CureTape Art has the same composition and adhesive strength as CureTape Classic and is also water resistant. CureTape Art is the only kinesiology tape manufactured using a unique printing process, in which the dye is completely incorporated into the fibres and does not lie on them. This keeps the cotton smooth, elastic and comfortable to the touch.



**Art Tattoo** 8717624162732  
**Art Tartan** 8717624162794  
**Art Tiger** 8717624163159  
**Art Zebra** 8717624163180  
**Art Leopard** 8717624163210  
**Art army Blue** 8717624164002






## CureTape® Beauty

**Materials:** Elastic cotton with a breathable, hypoallergenic acrylic adhesive.

**Dimensions:** Roll: 5 cm x 5 meters

**Displaybox:** 6 rolls

**Available colours:** 



**Beauty Beige**  
8717624165443

CureTape® Beauty is very suitable for use on facial skin. Face taping has become extremely popular to prevent or reduce the signs of aging by reducing wrinkles and lifting the skin. This is a wonderfully natural way to rejuvenate the facial appearance.




## MyCureTape® 3 rolls Art ON REQUEST

**Materials:** Elastic cotton with a breathable, hypoallergenic acrylic adhesive.

**Dimensions:** Roll: 5 cm x 2.5 meters

**Tube:** 3 rolls

**Available colours:** 

Each tube contains 2 rolls of CureTape Classic (5cm x 2.5m) and 1 roll of CureTape Art of 5cm x 2.5m, where you can choose from 4 trendy, contemporary prints. The various prints can be beautifully combined with the uni colors of CureTape Classic and are real eye-catchers! CureTape Art has the same composition as CureTape Classic, just with a nice added print.



**Tattoo Rebel**  
8717624165047



**Night Hunter**  
8717624165030



**Panther Power**  
8717624165016




**Safari Stripes**  
8717624165023

## CureTape® 2 rolls Beauty Tape

**Materials:** Elastic cotton with a breathable, hypoallergenic acrylic adhesive.

**Dimensions:** Roll: 5 cm x 5 meters

**Displaybox:** 2 rolls

**Available colours:** 



**Beauty 2 Roll Beige**  
8717624165801

CureTape® Beauty is very suitable for use on facial skin. Face taping has become extremely popular to prevent or reduce the signs of aging by reducing wrinkles and lifting the skin. This is a wonderfully natural way to rejuvenate the facial appearance. Kit includes x2 Facial silicone cups




## MyCureTape® 3 rolls Sports ON REQUEST

**Materials:** Elastic viscose tape with a breathable, hypoallergenic acrylic layer with 25% more adhesive strength than the Classic range.

**Dimensions:** Roll: 5 cm x 2.5 meters

**Tube:** 3 rolls

**Available colours:** 

Each tube contains 3 rolls of CureTape Sports (5cm x 2.5m) in various modern colours. The tubes are available in 3 variants. Extremely suitable during sports or in tape applications where stronger adhesion is required. The Sports range is made of elastic viscose and has an adhesive layer with 25% more adhesive strength than our cotton CureTape variants.

Packed in a practical, reusable tin tube! CureTape Sports is made of elastic viscose, it is water-resistant, has a hypoallergenic, acrylic-based adhesive layer and 25% more adhesive power than the regular CureTape classic variant. Viscose is just as comfortable to wear as cotton, and has a silkier and more supple feel and dries quickly.



8717624165467



8717624165764




8717624164965

## CureTape® Gentle

**Materials:** Elastic cotton with a breathable, hypoallergenic acrylic adhesive.

**Dimensions:** Roll: 5 cm x 5 meters

**Displaybox:** 6 rolls

**Available colours:** 




**Gentle Beige**  
8717624165689

CureTape® Gentle is very suitable for use on people with very sensitive skin, on young children and older people. CureTape® Gentle can also be used for tape applications in sensitive areas, such as on thin skin.

## CureTape® Self-Taping Book

The book is intended for people who want to tape an injury or complaint themselves. In this book, we have selected a number of common injuries for which kinesiology tape can help in case it is applied with the right technique. These taping applications can be used for sports or for carrying out your daily activities pain-free. Tape applications are supported by QR codes with a link to an instructional video



Author: Christina Peter

Language: English

More than 200 step-by-step instruction photos and illustrations

Year of release: 2022

Pages: 176

Self-Taping Guide  
9789081936330

## CureTape® Accessories



CureTape Clean Spray  
8717624161285



ON REQUEST

Dispenser  
8717624161001



Soft Touching  
Scissors  
8719172260717




ON REQUEST

High-Quality  
Scissors  
8719172260205

## MTC® Instruction Book

Written by MTC teacher Harry Pijnappel, this kinesiology taping method book has 192 pages with full-colour images. This manual is written in English especially for (para)medical professionals. The first chapter describes and substantiates the effects of the tape. Topics such as swelling, pain, itching and lymph are discussed. Subsequent chapters show practical applications for kinesiology taping. Basic taping techniques are discussed as well as numerous pathologies. Available in different languages.



Author: Harry Pijnappel

Language: English, Dutch, Spanish, German, French,  
Polish and Russian

Over 1300+ full colour images


Year of release: 2016

Pages: 192

MTC Handbook  
9789081397803

## MTC® Lymphatic Taping Book

This MTC® book on Treatment of Edemas, Hematomas and Scars by David Blow covers Decompressive Taping Techniques for Lymphatic and Vascular Rehabilitation and is a must-have if you are a health professional that applies kinesiology taping.



Author: David Blow

Language: English

Over 1300+ full colour images

Year of release: 2013

Pages: 413

Lymph taping Book  
9788870515992

# About CrossLinq®

CrossLinq® are non-elastic, grid-shaped tapes that are produced in our own factory in South Korea. CrossLinq has a skin-friendly acrylic adhesive layer. These tapes are mainly used on (local) pain points, trigger points and acupuncture points.

When applying to pain points you can easily apply the CrossLinq® yourself, however, application to trigger points and acupuncture points requires specific knowledge.



**CrossLinq Single Pack  
Small**  
(240 patches)  
8717624163012



**CrossLinq Single Pack  
Medium**  
(160 patches)  
8717624160998



**CrossLinq Single Pack  
Large**  
(40 patches)  
8717624163036

## CrossLinq® Patches


<b>Materials:</b>	Non-elastic, self-adhesive, waterproof, grid-shaped patches with a skin-friendly acrylic adhesive layer.
<b>Dimensions:</b>	Small: 1.4 x 1.8cm: 20 sheets x 12 crosslinqs Medium: 2.1 x 2.7cm: 20 sheets x 8 crosslinqs Large: 4.4 x 5.2cm: 20 sheets x 2 crosslinqs
<b>Displaybox:</b>	10 single packs
<b>Available Sizes:</b>	Small/Medium/Large

## CrossLinq® Book

The CrossLinq handbook contains more than 100 practical examples of treatment, 364 pages with clear explanations of the method, and high-quality photos of all treatments. The handbook also contains the latest insights and provides a total overview of the applications that are possible with CrossLinq, conveniently in alphabetical order of the conditions. CrossLinq is used on (local) pain points, trigger points, and acupuncture points. CrossLinq tapes are grid-shaped, electrostatically charged tapes that remove blockages, reduce or even remove pain and ailments, painlessly and without medication.

**Author:** Olaf Kandt  
**Language:** English  
**Instruction book with more than 100 applications**  
**Year of release:** 2016  
**Pages:** 364

**CrossLinq Book**  
9789082582611





# About FASCIQ®

FASCIQ® is a brand with high quality and easy to use fascia massage tools such as trigger point tools, IASTM tools and cupping sets. Recommended by leading R&D institutes.

FASCIQ® products are used to treat the fascia, to bring about a loosening of the fascia (fascia release) and to maintain a flexible fascial network (fascia treatment). Regular use of these products can prevent injuries.

Fascia mainly consist of connective tissue in different forms and structures. They are located throughout the body and are interconnected. In recent years, we have learned more and more about the function of the fascial network. Good results have been achieved with a treatment directed specifically at the fascia.

Among other things, the fascial network has an effect on:

- Strength
- Movement and mobility
- The immune system
- The balance of various internal systems (lymphatic, digestive, hormonal and metabolism)
- The posture

Various products are used to treat the fascia, to bring about a loosening of the fascia (fascia release) and to maintain a flexible fascial network (fascia treatment). Regular use of these products can prevent injuries.

## Endorsed by Dr. Robert Schleip

Fascia massage tools are an excellent method of releasing constrictions of the muscle. All FASCIQ® products are recommended by Dr Robert Schleip, director of the Fascia Research Group.



<https://toomacsolutions.co.nz/products/massage-tools>



# FASCIQ® IASTM Tools

The FASCIQ® IASTM Tools (“Instrument-Assisted Soft Tissue Mobilization Tools” or in short “Instrument Assisted Massage Tools”) consist of a number of separate tools with multiple shapes and edges to increase their employability. In addition, FASCIQ® also offers two complete IASTM toolsets with various shaped tools. All IASTM Tools are packed in a sturdy, beautiful storage bag.

The smaller tools are generally used for the treatment of smaller parts of the body, smaller strokes and reaching superficial tissue structures. The larger and heavier tools are used to reach deeper and are suitable for the treatment of larger parts of the body. These IASTM tools help prevent joint wear and strains complaints among therapists.



**Kiss**  
8719172260472



**Mustache**  
8719172260465



**Razor**  
8719172260489



**Drop**  
8719172260434



**Wheeler**  
8719172260441



**Wave**  
8719172260748



**Grip**  
8719172260427



**Artist**  
8719172260458

## FASCIQ® Small & Large Handlebar

**Large Handlebar**

Weight 0.539 kg  
Dimensions 44.6 × 1.9 × 0.9 cm

**Small Handlebar**

Weight 0.153 kg  
Dimensions 31 × 1 × 0.5 cm



**Large Handlebar**  
8719172260649

**Small Handlebar**  
8719172260632

# FASCIQ® Small Tool Set

- |    |                         |         |
|----|-------------------------|---------|
| 1) | 5 mm x 1 cm x 31 cm     | 153 gr  |
| 2) | 3 mm x 1,8 cm x 15,5 cm | 53 gr   |
| 3) | 3 mm x 1,8 cm x 12,6 cm | 39,5 gr |
| 4) | 3 mm x 3 cm x 11 cm     | 31,5 gr |
| 5) | 3 mm x 1,4 cm x 14,5 cm | 37,5 gr |



**Small Tool Set**  
8719172260410

# FASCIQ® Large Tool Set

- |    |                         |        |
|----|-------------------------|--------|
| 1) | 9 mm x 1,9 cm x 44,6 cm | 539 gr |
| 2) | 6 mm x 6 cm x 13,6 cm   | 207 gr |
| 3) | 6 mm x 4,2 cm x 20,7 cm | 263 gr |
| 4) | 6 mm x 3,2 cm x 21,3 cm | 216 gr |
| 5) | 5 mm x 2,3 cm x 15,2 cm | 85 gr  |



**Large Tool Set**  
8719172260403

# FASCIQ® Cupping Sets



The FASCIQ® silicone cupping sets (or silicone massage cups) are plastic silicone-based cups. 100% food grade silicone



Cupping set of 4  
8719172260373

Silicone Cupping Set of 4 different transparent Cellulite cups, packed in a convenient, FASCIQ® fabric storage bag. The diameters of these cuppings (cellulite cups) are 4.5 cm / 6 cm / 7.5 cm and 10 cm.



Small and Large Cup  
8719172260533

- Large Cup: This round cupping has a diameter of 7 cm and is 8 cm tall.
- Cup Small: This round cupping has a diameter of 5.5 cm and is 5.5 cm tall.



Facial cupping set  
8719172260540

The Mini facial cups are suitable for the treatment of eye contours and wrinkles above the lips. The Medium facial cups are used for neck, décolleté and other areas of the face.

- 2 x Medi cup : ø 3,8 cm x 8 cm
- 2 x Mini cup : ø 1,5 cm x 5 cm



Trigger point pods  
8719172260687

One piece design for treating different target areas  
2 suction level: plunger & inversion

- 2 x Small cups – 2.25" (57 mm)
- 2 x Large cups – 2.63" (67 mm)
- FASCIQ EVA carry case with zipper

# FASCIQ® Complete Cupping Set



This cupping collection set contains 14 cups in 10 different diameters.

All cups are made of high-quality silicone with a special grip addition. As a result, the cups offer extra grip during application. Ideal in combination with the free fascia cream included in the package! The set is suitable for both daily practice and home use.



Complete cupping set  
8719172261080

# FASCIQ® Accessories



FASCIQ Adhesive Spray  
8717624164033



Wax  
8719172260311



Cream  
8717624162640



2.5 cm (12 rolls)  
8719172261059

# FASCIQ® Sport Tape

**Materials:** Non-elastic, breathable, microporous cotton fabric with a latex-free, hypoallergenic and super-sticky, adhesive layer

**Dimensions:** 3 available sizes

**Available Colour:** White

Sport Tape is used by professionals of competitive sports to gain extra protection and prevent injuries.

**Features:**

- High tensile strength fabric provides enough support
- Serrated (zigzag) edge and hand tearable.
- No skin irritation, no residues on skin when peel off
- Stabilize joints during vigorous physical activity



3.8 cm (8 rolls)  
8719172260755

5 cm (6 rolls)  
8719172261066





Toomac Solutions  
32D Poland Road, Glenfield, Auckland 0627  
Freephone: 0508 443 5347  
Email: [Enquiries@toomac.co.nz](mailto:Enquiries@toomac.co.nz)  
Website: [www.toomac.co.nz](http://www.toomac.co.nz)